


Xmaru1717

Service Manual

Leading the transition to Digital X-ray





© Copyright 2010, Value Added Technologies Co. Ltd.

All pages of this document contain proprietary and confidential information of Value Added Technologies Corporation and are intended for exclusive use by Value Added Technologies Corporation personnel or customers. Copying, disclosure to others or other use is prohibited without the express written authorization from Value Added Technologies Corporation. Please report any violations of requirement to Value Added Technologies Corporation.

Attention

For improvement of product performance, supplementation, or follow-up of information; the contents of this manual are subject to change without separate prior notice.

Please note that our company has neither responsibility for any accidents nor obligation to do free repair service for any damage of the equipment due to user's mistake, which resulted from failure to follow the contents in this manual. Make sure to be familiar with the safety precautions and usage procedures. Also note that the product may slightly differ from the contents of this manual depending on specification.

The following marks are used for the effective use of the product in this manual.



Attention, consult accompanying documents.



This is used to emphasize essential information. Be sure to read this information to avoid incorrect operation.



This indicates hazardous situation which, if not heeded, may result in minor or moderate injury to you or others, or may result in machine damage.



This indicates a potentially hazardous situation which, if not heeded, could result in death or serious injury to you or others.

Federal Law restricts this device to sale by or the order of a radiologist or any other practitioners licensed by the law of the state in which that person practices to use or order the use of the device.

Contents

Chapter 1	Introduction 1	
	Overview	1
	Product features	1
	Product specification	
Chapter 2	IP Address Set-up	3
	IP switching	
	Mac address set-up	
Chapter 3	Firmware Set-up	11
	FTP set-up	11

Introduction

Overview

Xmaru1717 is an X-ray image acquisition device that is based on flat-panel. This device should be integrated with an operating PC and an X-ray generator. It can do to utilize as digitalizing X-ray images and transfer for radiography diagnostic.

Product features



Figure 1: Xmaru1717

The surface of housing is made of the metal and X-ray window area is covered with carbon. Inner circuit of detector uses lead to protect it from X-ray and has a structure protecting the inner circuit from X-ray dose under the limits

- Based on a-Si TFT active matrix
- Wide image
- 14-bit digital output
- Easy integration

Product specification

Parameter	Unit	Description
Application	-	General radiography
Pixel Size	μm	143
Active Pixel Array	pixels	3072 x 3072
Active area	mm	439.3 x 439.3
Spatial Resolution* ¹	lp/mm	3
Dynamic Range	dB	> 73
Energy Range	kVp	40 ~ 150
A/D	bit	14 bit, 16834 gray scale
Dimension (W*H*D)	mm	500 x 497 x 45
Weight	kg	13.4
Data Output	-	100 Mbps
Power	V	AC100-120V/200-240V (50 / 60Hz)
Environmental Condition		Operation: Temperature: +10 ~ +40 °C Humidity: 20 ~ 75 % (no condensation) Press 70 ~ 106kPa Storage and transportation: Temperature: -25~ +55 °C Humidity: 10 ~ 95 % (no condensation) Press 70 ~ 106kPa

^{*1:} line pair resolution at over 10% MTF

Table 1: Product specification

IP Address Set-up

IP switching

Step 1




Figure 2 : Telnet connection

- Connect the detector with PC (Installation should be done)
- Click "START" and then "RUN"
- Input "telnet 192.168.1.80"

Step 2




Figure 3: Telnet login

Input

ID: "root" and

Password: "depss"

in Pop-up "telnet 192.168.1.80"

VATECH login: root

Password : depss (Password does not appear in window)

```
Telnet 192.168.1.80

VATECH login: root

Password:

[root@VATECH root]$vi /etc/config/ifcfg-eth0_
```

Figure 4: Input massage

- Detector is connected, when the message "[root@VATECH root]\$" appears
- Input the message
 - "vi/etc/config/ifcfg-eth0", adjacent to "[root@VATECH root]\$"

Step 4

```
Telnet 192.168.1.80

DEUICE=eth0
IPADDR=192.168.1.80
BROADCAST=192.168.1.255
NETMASK=255.255.255.0
GATEWAY=192.168.1.1
```

Figure 5 : IP address edit 1

• Move curser at under "0" of "80"

```
Telnet 192.168.1.80

DEUICE=eth0
IPADDR=192.168.1.80_
BROADCAST=192.168.1.255
NETMASK=255.255.255.0
GATEWAY=192.168.1.1
```

Figure 6: IP address edit 2

• Input "a" and then "80_" appears. And then " - Insert - " appears.

Step 6

```
Telnet 192.168.1.80

DEUICE=eth0
IPADDR=192.168.1.8_
BROADCAST=192.168.1.255
NETMASK=255.255.255.0
GATEWAY=192.168.1.1
```

Figure 7: IP address edit 3

Push the "Delete" button, and then "8_" appear.

```
Telnet 192.168.1.80

DEUICE=ethØ
IPADDR=192.168.1.81
BROADCAST=192.168.1.255
NEIMASK=255.255.255.0
GATEWAY=192.168.1.1
```

Figure 8: IP address edit 4

Input "1"
When inputting another IP address, it's possible to input another IP address using the "delete" button and inputting the number directly. And to launch dual detectors is possible, the number is variable among "0~255" but not be repeated.

Step8




Figure 9: IP address edit 5

Push "ESC" button and input ":".
 Push "ESC" button and then the message " – Insert –" disappears, and push ":" and then ":" appears..

```
Telnet 192.168.1.80

DEUICE=eth0
IPADDR=192.168.1.81
BROADCAST=192.168.1.255
NETMASK=255.255.25.0
GATEWAY=192.168.1.1
```

Figure 10: IP address edit finish

• Push "wq!" and push "Enter" button and then IP address is set. If you don't want to set IP address, input "q!" and repeat again from "step 3".

Mac address set-up

After finishing the IP address setting, to distinguish Ethernet Mac address should be changed.

Step 1

```
Telnet 192.168.1.80

#! /bin/sh

# /etc/rc$.d/10network

# Bring up the network and interfaces.

[ -f /etc/config/network ] && . /etc/config/network

# Setting the hostname
echo "Setting hostname..."
hostname $HOSTNAME

# Setting forwarding and defrag
if [ -d /proc/sys/net/ipv4 ]; then
    if [ "$FORWARD_IPV4" = "no" ]; then
    value=0
    message="disabled"
    else
    value=1
    message="enabled"
    fi echo "$value" > /proc/sys/net/ipv4/ip_forward
    echo "IPv4 packet forwarding is $(message)..."
"/etc/rc.d/init.d/network" line 1 of 89 --1%---
```

Figure 11: Mac address edits 1

 Adjacent to the message "[root@VATECH root]\$", input the message "vi /etc/rc.d/init.d/network" and then it will appear as below.

Step 2

Figure 12 : Mac address edits 2

- Move curser to the message
 - "# Bringing up the local network"

Figure 13: Mac address edits 3

- Input "I", then the message " Insert " appears and input will be possible.
- Input the message "ifconfig eth0 hw ether 00:64:69:73:6x".
 The value of X is available among "1~E"

EX) First Detector → 00:64:67:69:73:61

Second Detector → 00:64:67:69:73:62

Step 4

Figure 14: Mac address edits 4

Figure 15: Mac address edits 5

- After inputting value is finished, set again with the steps from 3.a ~ 3.b
- Push "ESC" button, input ": ".
- After input "wq!", push "Enter" button and then IP address is set. If setting is wrong, input "q!" and redo from the step 1.

Firmware Set-up

FTP set-up

Step 1: Free FTP download

http://filezilla-project.org/


http://sourceforge.net/projects/filezilla






Figure 16: Free FTP downloads URL

Run "FileZilla_Server-N_N_NN.exe" file








Figure 17: FTP installation 1




Figure 18:FTP installation 2

Step 2: FTP set-up step by step

A. Setting

[Menu Bar] \rightarrow [Edit] \rightarrow [Settings]




Figure 19: FTP set up

General settings




Figure 20: General settings

Passive mode settings




Figure 21: Passive mode settings

B. Users

[Menu Bar] \rightarrow [Edit] \rightarrow [Users]




Figure 22: Users settings

General → Add






Figure 23 : User Add

Shared folders → Add






Figure 24: Shared folders Add

Step 3: Detector connection





Figure 25: Detector connection check

Step 4: Telnet connection



```
Telnet 192.168.1.80
 root@VATECH home]$
root@VATECH home]$
root@VATECH home]$pwd
 /home
  root@VATECH home]$
  root@VATECH home]$1s -a1
                    2 root
15 root
                                                               O Dec
                                                                              2030
 irwxr-xr-x
                                                                         4
                                       root
                                                               0 Jan
                                                                               1970
 drwxr-xr-x
                                       root
                                                                               1970 .cfg
                                       root
                                                              84 Jan
                         root
                                                               4 Dec
                         root
                                                                               2030 .el_fld
                                       root
                                                                         4
                         root
                                       root
                                                               4 Dec
                                                                              2030 .el_org
                                                       338199 Dec
                                                                              2030 davinci
  rwxr-xr-x
                      1 root
                                       root
                                                                         4
rw-r--r-
root@VATECH home]$

root@VATECH home]$ps

PID Uid VmSize Stat Command

1 root 1072 S init

2 root S [kevento

3 root S [kswapd

4 root S [kswapd

5 root S [bdflus

5 root S [kupdat

8 [mtdblo
                                                              41 Dec 4 2030 version.txt
  rw-r--r-
                      1 root
                                       root
                                           [keventd]
                                           (ksoftirqd_CPUO)
                                           [kswapd]
                                           bdf lush]
                                           [kupdated]
                                           [mtdblockd]
                                         [jffs2_gcd_mtd3]
devfsd /dev
      8 root
                        1120 S
1840 S
1340 S
1116 S
1732 S
1056 S
1732 S
21972 S
21972 S
21972 S
21972 S
21972 S
     22 root
39 root
                                          /bin/bash /etc/rc.d/rc 5
     66 root
                                          inetd
     73 root
74 root
77 root
78 root
79 root
80 root
                                          initlog -c /etc/rc.d/rc5.d/S991ocal start
                                          /bin/sh /etc/rc.d/rc5.d/S991ocal start
/bin/utelnetd -d
                                          /bin/sh /usr/bin/run
                                          /home/davinci
                                          /home/davinci
     81 root
                                          /home/davinci
                                         /home/davinci
         root
     85 root
                                          -sh
                           1364 R
     95 root
                                          ps
  root@VATECH home]$
   root@VATECH home]$_
```

Figure 26: Telnet connection

- "[root@VATECH home ~]# telnet 192.168.1.80" (connect detector)
- Vatech login : root
- Password : depss
- "[root@VATECH home ~]#cd home" (Move to home directory)
- "[root@VATECH home ~]#pwd" (Check current directory)
- "[root@VATECH root]\$ps" (Check current process information)

Step 5: Firmware download

```
Telnet 192.168.1.80
[root@VATECH home]$
[root@VATECH home]$1s
version.txt
[root@VATECH home]$
Proot@UATECH home l$ftp 192.168.1.100
Connected to 192.168.1.100.
220-FileZilla Server version 0.9.29 beta
220-written by Tim Kosse (Tim.Kosse@gmx.de)
220 Please visit http://sourceforge.net/projects/filezilla/
Name (192.168.1.100:root): root
331 Password required for root
Password:
230 Logged on
Remote system type is UNIX.
ftp> ?Invalid command
ftp> ls
200 Port command successful
150 Opening data channel for directory list.
                               1645568 Oct 07 2008 davinci_080617(2).tar
157090 Sep 29 2008 davinci_080617.zip
 rw-r--r-- 1 ftp ftp
rw-r--r-- 1 ftp ftp
 -rw-r--r-- 1 ftp ftp 338199 Dec 05 17:03 davinci_D007-081205
226 Transfer OK
ftp> ftp> bin davinci_D007-081205
200 Type set to I
ftp> ftp> get davinci_D007-081205
200 Port command successful
150 Opening data channel for file transfer.
226 Transfer OK
338199 bytes received in 1.6 seconds (211465 bytes/s)
ftp> ftp>
```

Figure 27: Firmware download

- "[root@VATECH home ~]\$ ftp 192.168.1.100" (connect FTP)
- "[root@VATECH home ~]\$ ftp user ID (Ex:root), PW(Ex:1234) (davinci_D007-081205, firmware will be uploaded)
- "ftp> bin davinci_D007-081205" (convert "davinci_D007-081205" to binary file)
- "ftp> get davinci_D007-081205" (download "davinci_D007-081205" to detector)
- "ftp> bye" (end ftp)

Step 6: Firmware set-up

```
Telnet 192.168.1.80
[root@VATECH home]$
[root@VATECH home]$1s
davinci_D007-081205 version.txt
[root@VATECH home]$
[root@VATECH home]$chmod 755 davinci_D007-081205
[root@VATECH home]$
[root@VATECH home]$1s -al
                                            Ø Dec 4 02:54
drwxr-xr-x
              2 root
                           root
             15 root
                                           Ø Jan 1
84 Jan 1
                                                       1970
drwxr-xr-x
                           root
              1 root
                                                       1970 .cfg
 rw-r----
                           root
                                             4 Dec
                                                       2030 .el_fld
               1 root
                           root
               1 root
                           root
                                             4 Dec
                                                       2030 .el_org
 rw-r
                                       338199 Dec
                                                       02:54 davinci_D007-081205
 PWXP-XP-X
               1 root
                           root
                                                    4
                                            41 Dec
                                                    4 2030 version.txt
               1 root
                           root
[root@VATECH home]$
[root@VATECH home]$
```

Figure 28: Firmware set-up

"[root@VATECH home ~]#chmod 755 davinci_D007-081205"

```
(set authority: User :7 (read, write, execute)

Group :5 (read, execute)

Other 5 (read, execute)
```

```
Telnet 192.168.1.80

davinci_D007-081205 version.txt

[root@VATECH home]$
[root@VATECH home]$mv davinci_D007-081205 davinci
[root@VATECH home]$
[root@VATECH home]$
[root@VATECH home]$15 -1
-rwxr-xr-x 1 root root 338199 Dec 4 1916 davinci
-rw-r-r-- 1 root root 41 Dec 4 2030 version.txt
[root@VATECH home]$
[root@VATECH home]$
```

Figure 299: Convert file name

"[root@VATECH home ~]# mv davinci_D007-081215 davinci"

(convert file name "Davinci_007-081215" to davinci)



```
Telnet 192.168.1.80
VATECH login: root
Password:
[root@UATECH root]$ps
 PID Uid
               UmSize Stat Command
   1 root
                 1072 S
                            init
    2 root
                            [keventd]
                            [ksoftirqd_CPU0]
   3 root
    4 root
                            [kswapd]
   5 root
                            [bdflush]
                           [kupdated]
   6 root
   7 root
                            [mtdblockd]
   8 root
                           [jffs2_gcd_mtd3]
   22 root
                 1120 S
                           deufsd /dev
   39 root
                 1840 S
                           /bin/bash /etc/rc.d/rc 5
   66 root
                 1340 S
                            inetd
   73 root
                 1116 S
                           initlog -c /etc/rc.d/rc5.d/$991ocal start
   74 root
                 1732 S
                           /bin/sh /etc/rc.d/rc5.d/S99local start
   77 root
                 1056 S
                           /bin/utelnetd -d
                            /hin/sh/usr/bin/run
   78 root
                 1732 S
   79 root
                21972 S
                           /home/davinci
   80 root
                21972 S
                           /home/davinci
   81 root
                21972 S
                           /home/davinci
   82 root
                21972 S
                           /home/davinci
   84 root
                 1772
   87 root
                 1364 R
                           ps
[root@VATECH root]$
[root@VATECH root]$
```

Figure 30: Check "davinci" process updated

Check "davinci" process updated

Step 7 : Firmware set-up




Figure 31: Check firmware version



Service Manual

Leading the transition to Digital X-ray





VATECH Co., Ltd.